



Online store [WWW.LYNXGEAR.LV](http://WWW.LYNXGEAR.LV)  
Company: SIA LATLYNX  
info@lynxgear.lv, +371 26159668  
Reg. No. 45403038162, VAT Reg. No.  
LV45403038162,  
Legal address: Baznicas iela 8,  
Aizkraukle, Aizkraukles novads,  
LV-5101, LATVIJA  
Bank Account: LV87UNLA0050021387494  
AS SEB Bank, SWIFT: UNLALV2X

---

## ORDER CONFIRMATION

---

Not Rated Yet

**Sales price 690,00 €**

Salesprice with discount

Incl. VAT 21%: 119,75 €



### Description

A tried and trusted design, used as a staple item by Emergency Services and Rescue Organisations for many years. Its main function, apart from transporting the injured to safety, is to do so in such a way that protects the spine where an accident is suspected of causing damage. Instead of a flat, plain bed, the stretcher incorporates a frame of hinged panels. In their unlocked position, each panel is capable of being swivelled, (rather like a cellar trapdoor). Without disturbance to the patient's position, the stretcher can be manoeuvred in such a way that allows each of these panels or flaps to slide underneath the body. Once in its final position, all the plates are locked into place and the patient can be safely moved. Instead of being lifted onto a traditional stretcher. One vital feature is the way the stretcher maximises patient immobilisation during transit. Such details can actually make the difference between recovery and permanent spinal damage. Another important feature of the Relequip 2-Piece Rescue Stretcher is that when locked in place, the panels that form the bed leave a narrow



Online store [WWW.LYNXGEAR.LV](http://WWW.LYNXGEAR.LV)  
Company: SIA LATLYNX  
info@lynxgear.lv, +371 26159668  
Reg. No. 45403038162, VAT Reg. No.  
LV45403038162,  
Legal address: Baznicas iela 8,  
Aizkraukle, Aizkraukles novads,  
LV-5101, LATVIJA  
Bank Account: LV87UNLA0050021387494  
AS SEB Bank, SWIFT: UNLALV2X

---

## **ORDER CONFIRMATION**

---

gap that runs the full length of the stretcher. This allows Emergency Staff to carry out Spinal X-Ray without moving the patient.

For gentle patient manoeuvring without rolling or lifting

Open center allows patient to be x-rayed

Lightweight aluminium

Adjustable length

Easy to clean

Open 6cmH x 44cmW x 190cmD

Closed 6cmH x 44cmW x 120cmD

Great Britain.